



Intended Use

FoundationOne® Liquid CDx is a next generation sequencing based in vitro diagnostic device that analyzes 324 genes. Substitutions and insertion and deletion alterations (indels) are reported in 311 genes, copy number alterations (CNAs) are reported in 310 genes, and gene rearrangements are reported in 324 genes. The test also detects tumor fraction and the genomic signatures blood tumor mutational burden (bTMB) and microsatellite instability high (MSI-H) status. FoundationOne® Liquid CDx utilizes circulating cell-free DNA (cfDNA) isolated from plasma derived from the anti-coagulated peripheral whole blood of cancer patients. The test is intended to be used as a companion diagnostic to identify patients who may benefit from treatment with targeted therapies in accordance with the approved therapeutic product labeling. Additionally, FoundationOne® Liquid CDx is intended to provide tumor mutation profiling to be used by qualified health care professionals in accordance with professional guidelines in oncology for patients with malignant neoplasms.

A negative result from a plasma specimen does not mean that the patient's tumor is negative for genomic findings. Patients who are negative for genomic findings should be reflexed to routine biopsy and their tumor mutation status confirmed using an FDA-approved tumor tissue test, if feasible.



Summary of Analytical Sensitivity and Specificity¹

Results from our Limit of Detection (LoD) study are shown below, indicating the median variant allele frequency, tumor fraction or unstable loci at which the test has shown 95% probability of detection. Please refer to our product labeling for a list of the 75 genes baited for enhanced sensitivity and complete product specifications.

ALTERATION TYPE	BAIT SET REGION	MEDIAN LIMIT OF DETECTION (LOD)
Short Variants	Enhanced Sensitivity	0.40% VAF
	Standard Sensitivity	0.82% VAF
Rearrangements	Enhanced Sensitivity	0.37% VAF
	Standard Sensitivity	0.90% VAF
Copy Number Amplification	NA	21.7% TF
Copy Number Loss	NA	30.4% TF
MSI	NA	0.8% Unstable loci
bTMB (component indels)	NA	1.00% VAF
bTMB (component subs)	NA	1.00% VAF

VAF = variant allele frequency; TF = tumor fraction
The accuracy of %VAF / %TF have not been analytically validated

In our Limit of Blank study, which evaluated variant calling in healthy donors, 1,735 unique variants were included in the analysis for a total of 137,065 data points. A total of 18 false positives were observed across 4 unique short variants. The LoB was determined to be the ideal value of zero for short variants, rearrangements and CNAs. The false positive rate was shown to be 0% for rearrangements and CNAs and 0.013% (-1 in 8,000) for short variants (substitutions and indels).



FoundationOne® Liquid CDx Gene List¹

As a professional service, FoundationOne® Liquid CDx interrogates 324 genes, including 309 genes with complete exonic (coding) coverage and 15 genes with only select non-coding coverage (indicated with an*); **75 genes (indicated in bold) are captured with increased sensitivity** and have complete exonic (coding) coverage unless otherwise noted. The test also detects tumor fraction and the genomic signatures blood mutational burden (bTMB) and microsatellite instability high (MSI-H) status.

ABL1 [Exons 4-9]	ALOX12B	ASXL1	BAP1	BCR* [Introns 8, 13, 14]	BRIP1	CASP8
ACVR1B	AMER1 (FAM123B)	ATM	BARD1	BRAF [Exons 11-18, Introns 7-10]	BTG1	CBFB
AKT1 [Exon 3]	APC	ATR	BCL2	BRCA1 [Introns 2, 7, 8, 12, 16, 19, 20]	BTG2	CBL
AKT2	AR	ATRAX	BCL2L1	BRCA2 [Intron 2]	BTK [Exons 2, 15]	CCND1
AKT3	ARAF [Exons 4, 5, 7, 11, 13, 15, 16]	AURKA	BCL2L2	BRD4	C11orf30 (EMSY)	CCND2
ALK [Exons 20-29 Introns 18,19]	ARFRP1	AURKB	BCL6		C17orf39 (GID4)	CCND3
	ARID1A	AXINI	BCOR		CALR	CCNE1
		AXL	BCORL1		CARD11	CD22

(FoundationOne® Liquid CDx Gene List continued)

CD70	ERBB2	FOXL2	KLHL6	NF1	PPARG	SMAD2
CD74* [Introns 6-8]	ERBB3 [Exons 3, 6, 7, 8, 10, 12, 20, 21, 23, 24, 25]	FUBP1	KMT2A (MLL) [Introns 6, 8-11, Intron 7]	NF2	PPP2R1A	SMAD4
CD79A	ERBB4	GABRA6	NFE2L2	NFKBIA	PPP2R2A	SMARCA4
CD79B	ERCC4	GATA3	NFKB1	PRDM1	PRKAR1A	SMARCB1
CD274 (PD-L1)	ERG	GATA4	NKX2-1	PRKAR1A		SMO
CDC73	ERRF1	GATA6	NOTCH1	PRKCI		SNCAIP
CDH1	ESR1 [Exons 4-8]	GNA11 [Exons 4, 5]	LTK	PTCH1		SOCS1
CDK12	ESR1 [Exons 4-8]	GNA13	LYN	PTEN		SOX2
CDK4	ETV4* [Intron 8]	GNAQ [Exons 4, 5]	MAF	NOTCH3	PTPN11	SOX9
CDK6	ETV5* [Introns 6,7]	GNAS [Exons 1, 8]	MAP2K1 (MEK1) [Exons 2, 3]	NPM1 [Exons 4-6, 8, 10]	PTPRO	SPEN
CDK8	ETV6* [Introns 5,6]	GRM3	MAP2K2 (MEK2) [Exons 2-4, 6, 7]	NRAS [Exons 2, 3]	QKI	SPOP
CDKN1A	EWSR1* [Introns 7-13]	GSK3B	MAP2K4	NSD3 (WHSC1L1)	RAC1	SRC
CDKN1B	EZH2 [Exons 4, 16, 17, 18]	H3F3A	MAP3K1	NT5C2	RAD21	STAG2
CDKN2A	EZR* [Introns 9-11]	HDAC1	MAP3K13	NTRK1 [Exons 14, 15, Introns 8-11]	RAD51	STAT3
CDKN2B	FAM46C	HGF	MAPK1	NTRK2 [Intron 12]	RAD51B	STK11
CDKN2C	FANCA	HNFA	MCL1	NTRK3 [Exons 16, 17]	RAD51C	SUFU
CEBPA	FANCC	HRAS [Exons 2, 3]	MDM1	NUTM1* [Intron 1]	RAD51D	SYK
CHEK1	FANCG	HSD3B1	MDM4	P2RY8	RAD52	TBX3
CHEK2	FANCL	ID3	MED12	PALB2	RAD54L	TEK
CIC	FANCL	IDH1 [Exon 4]	MEF2B	PARK2	RAF1 [Exons 3, 4, 6, 7, 10, 14, 15, 17, Introns 4-8]	TERC* {ncRNA}
CREBBP	FAS	IDH2 [Exon 4]	MEN1	PARP1	RARA [Intron 2]	TERT* {Promoter}
CRKL	FBXW7	IGF1R	MERTK	PARP2		TET2
CSF1R	FGF10	IKBKE	MET	PARP3	RB1	TGFBR2
CSF3R	FGF12	IKZF1	MITF	PAX5	RBM10	TIPARP
CTCF	FGF14	INPP4B	MKMK1	PBRM1	REL	TMPRSS2* [Introns 1-3]
CTNNA1	FGF19	IRF2	MLH1	PDCD1 (PD-1)	RET [Introns 7, 8, Exons 11, 13-16, Introns 9-11]	TNFAIP3
CTNNB1 [Exon 3]	FGF23	IRF4	MPL [Exon 10]	PDCD1LG2 (PD-L2)		TNFRSF14
CUL3	FGF3	IRS2	MRE11A	PDGFRA [Exons 12, 18, Introns 7, 9, 11]		TP53
CUL4A	FGF4	IRS4	MSH2 [Intron 5]	PDGFRB [Exons 12-21, 23]		TSC1
CXCR4	FGF6	JAK1	MSH3	PDK1		TSC2
CYP17A1	FGFR1 [Introns 1, 5, Intron 17]	JAK2 [Exon 14]	MSH6	PIK3C2B		TYRO3
DAXX	FGFR2 [Intron 1, Intron 17]	JAK3 [Exons 5, 11, 12, 13, 15, 16]	MST1R	PIK3C2G		U2AF1
DDR1	FGFR3 [Exons 7, 9 (alternative designation exon 10), 14, 18, Intron 17]	JUN	MTAP	PIK3CA Exons 2, 3, 5-8, 10, 14, 19, 21 (Coding Exons 1, 2, 4-7, 9, 13, 18, 20)		VEGFA
DDR2 [Exons 5, 17, 18]	FGFR4	KDM5A	MTOR [Exons 19, 30, 39, 40, 43-45, 47, 48, 53, 56]	PIK3CB		VHL
DIS3	FH	KDM5C	MUTYH	PIK3R1		WHSC1
DNMT3A	FLCN	KDM6A	MYB* [Intron 14]	PIM1		WT1
DOT1L	FLT1	KDR	MYC [Intron 1]	PMS2		XPO1
EED	FLT3 [Exons 14, 15, 20]	KEAP1	MYCL (MYCL1)	POLD1		XRCC2
EGFR [Introns 7, 15, 24-27]		KEL	MYCN	POLE		ZNF217
EP300		KIT [Exons 8, 9, 11, 12, 13, 17, Intron 16]	MYD88 [Exon 4]			ZNF703
EPHA3			NBN			
EPHB1						
EPHB4						

1. Data on file M-GB-00001408



© 2020 Foundation Medicine, Inc. | Foundation Medicine® and FoundationOne® are registered trademarks of Foundation Medicine, Inc.

Roche Products Limited
Hexagon Place, 6 Falcon Way, Shire Park, Welwyn Garden City, Hertfordshire AL7 1TW.
Registered in England 100674
medinfo.uk@roche.com 0800 3281629
Roche Products Limited is the sole licensed distributor of Foundation Medicine, Inc. products in the UK

M-GB-00004362 Date of Preparation July 2021 RAL-0106-01